

RPV-1 Generation 3 Polymer Viscometer

Application: Determination of the reduced viscosity (also known as viscosity number) and *K*-value of PVC resins dissolved in **Cyclohexanone in accordance with ISO 1628 Part 2**



RPV-1 (2) Gen. 3 RSS

Supplied with Auto sampler, optional sample preparation (iSP-1) and reaction block (SD-1BM not shown)

RPV-1 (2) Gen. 3 Auto

Supplied with Auto sampler, sample preparation (iSP-1) and single (or optionally dual) reaction block (SD-1BM not shown)

